

Monocular Microscope



Technical Specification:

MAIN BODY	Ergonomically designed inclined body, all die casted sturdy body
VIEWING HEAD	45 degree inclined monocular head
OBJECTIVES	Anti-fungal, antireflection coated 4x, 10x, 40x and 100x oil immersion aemi-plan achro din objectived, fully brass made with clear lenses and nicely designed external look
LIGHT SOURCE	high quality L.E.D. illumination & SMPS based 100V 220V wide voltage range
EYE-PIECES	Anti-futi reflection coated - 10x wide-field and 15x wide-field pairs
STAGE	Graduated mechanical rectangular stage.
NOSEPIECE	Outward quadruple brass made nosepiece with stop clicks.
MAGNIFICATION SETTINGS	8 magnification settings: 40X, 60X, 100X, 150X, 400X, 600X, 1000X & 1500X
CONDENSER	ABBE 1.25 NA IRIS diaphragm condenser with aluminium made IRIS diaphragm & swing in-swing out IRIS diaphragm
PACKING AND ACCESSORIES	In box content-main body, binocular head, 4 objectives, 4 eye-pieces, power cable, dust cover, cleaning cloth, 50 blank slides and cover slips. Packed in styroform box and cardboard box with instruction manual.

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.