

NSG-3000FT

Nitrogen Spray GUN



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Nitrogen Spray Guns

N2 Spray Gun Assembly - Straight Line

TEQCOM's PTFE Nitrogen Guns deliver the optimum in anti-corrosion protection and long lasting duty whenever harsh chemicals are used in the vicinity of N2 dispensing/drying. Manufactured of PTFE and poly (like the companion piece ... DI Water Spray Guns) they resist acid fume attack. Moreover, the N2 guns' unique design eliminates costly nitrogen leaks



Item#	Description	Packed
TA-N2-3100	N2 Spray Gun, Straight Assembly	One per carton

(which occur in conventional guns around the depress button). N2 guns are available individually or as assemblies with a hose. This Standard N2 Spray Gun is supplied with a 7 foot hose and fittings for N2 supply hookup, filter housing and one TACODISC-10 Filter.

►► N2 Spray Gun Assembly

N2 Gun in spray gun configuration, with 7 foot coiled hose and fittings for N2 supply hookup, filter housing and one each disposable TACODISC Filter and front or rear trigger. PTFE NITROGEN GUNS deliver the optimum in anti-corrosion protection and long lasting duty wherever harsh chemicals are used in the vicinity of N2 dispensing/drying. Manufactured of PTFE and poly (like the companion piece . . . our DI Water Spray Guns), they resist acid fume attack. Moreover, our N2 guns' unique design eliminates costly nitrogen leaks (these leaks usually occur in conventional guns around the depress button). N2 Guns are available individually or as assemblies with our popular coiled hose.

Item#	Description	Packed
TA-N2-3200FT	N2 Spray Gun Assembly, Front Trigger	One per carton
TA-N2-3200RT	N2 Spray Gun Assembly, Rear Trigger	One per carton



▶▶ Nitrogen Spray Guns Replacement Parts for N2 Spray Guns

Replacement Filter 0.45 Micron for Nitrogen Spray Gun

Item#	Description	Packed
TA-CODISC-10	0.45 Micron Filter for N2 Spray Gun	10 per carton



OD Hose - 1/2" for use with TA-N2-2000 Spray Gun

Item#	Description	Packed
TA-17001	OD 1/2" Hose	One per carton

Nitrogen Spray Gun Hose w/Fittings

Item#	Description	Packed
TA-18003	Nitrogen Spray Gun Hose 1/4" Fittings on both ends. 1/8" ID, 1/4" OD	One per carton




Nitrogen Spray Gun Hose

Item#	Description	Packed
TA-18004	Nitrogen Spray Gun Hose 7 foot flexible coiled hose. 1/2" ID, 1/4" OD	One per carton

Spray Gun Illustration

N2 Spray Gun Assembly



Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide