

OSI - 3112

Shaking Incubator



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Features:

- Audio and visual alarm.
- Power off recovery function.
- Soft shaking mode, stable and reliable operation.
- Efficient drive motor and cooling controller, low energy consumption.
- Transparent window and built-in lighting design convenient for observation samples.
- Accessories: UV lamp(standard accessories); Humidity controller, Inlet air pipe, Light system (optional accessories).

►► Specifications:

Model	OSI-3112	OSI-3212
Controller:	PID	
Display:	LCD	
Circulation Mode:	Forced convection	
Shaking Mode:	Rotary shaking	
Drive Mode:	Unishaft drive	
Rotating Speed:	30~300rpm	
Rotating Accuracy:	±1rpm	
Vibration Range:	Φ25mm	
Temp. Range:	RT+5~60°C	4~ 60°C
Temp. Accuracy:	±0.1°C	
Temp. Uniformity:	±1°C	
Timing Range:	0~999.59h	
Shaking Plate Size:	500*350mm	
Shaking Plate Quantity:	2 pcs	
Max. Capacity:	1000ml*12 or 500ml*24 or 250ml*30 or 100ml*48 or 50ml*48	
Standard Capacity:	500ml*12 ,250ml*15	
Internal Material:	Stainless steel	
External Material:	Cold-rolled steel with anti-bacteria power coating	
Refrigerant:	/	CFC-Free 134A
Cooling Consumption:	/	200W
Heating Consumption:	500W	
Safety Function	Audio and visual alarm for over-temperature;	
	Alarm for over-speeds;	
	Independent over-temperature protector;	
	Shutdown protection for door opening;	
	Leakage or over-current protection.	
Consumption:	600W	830W
Power Supply:	AC110V/220V±10%,50/60HZ	
Internal Size(W*D*H)(mm)	600*440*650	
External Size(W*D*H)(mm)	700*600*1250	
Package Size(W*D*H)(mm)	800*700*1400	
Gross Weight(kg)	232	

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide