

PC - 3020

Plate Centrifuge



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Plate Centrifuge is a specially designed personal-sized centrifuge for multi-plates, such as PCR plate, ELISA plate, cell culture plate and the deep well plate.

PC-3020 can carry up to 4xELISA plates or 2xPCR plates or 2xdeep well plates at the same time. The max. speed is 2000rpm with the centrifugal force(RCF) of 600xg.

PC-3020 runs quietly, driven directly by the high force torque motor. It shows time and speed intuitively in digits. The loading of the samples is easy and convenient. It has a fast ramping rate of speed and has the “Spin” instantaneous centrifugal function.

PC-3020 is your excellent assistant in the lab as a personal plate centrifuge with the minimum size but the strong capability.

►► Features :

- Strong force- directly driven by torque motors;
- Max. speed up to 2000rpm, Max. RCF up to 600xg;
- Setting time and speed by pressing buttons easily;
- “Spin”- the instantaneous centrifugal function;
- Microcomputer controlled unit, time and speed digital display;
- Safety locker, automatically stop while opening the lid;
- Small footprint and light weight.



►► Specification :

Model	PC-3020
Max. speed	2000 rpm
Max. RCF	600g
Sample capacity	4xELISA plates or 2xPCR plates or 2xDeep well plates
Speed setting range	500-2000 rpm
Noise	<60dB(A)
Time setting range	1-99 mins, continuous
Max. heating tempt.	RT+12°C
Control system	Microcomputer
Driving system	Directly driven by high force torque motor
Dimension	31cm x 38cm x 30cm
Power Supply	AC 100V~110V/200V~220V, 50/60Hz, 150W
Net weight	12Kg

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide