

Penta Head Microscope



Main Features:

PENTA HEAD MICROSCOPE Incorporates in innovative optical and structural design idea, facilitating group discussions and consultation by five observers during important meetings, conferences and seminars. Penta Head is equipped with five observation head, inclined at 45°/30° for comfortable viewing, including the central Trinocular Head to accept optional Photomicrography and Video Equipments. The built-in LED Green Light Pointer, with adjustable brightness can be moved in the field of view during group discussion. Supplied as above with following infinity corrected plan optics.

Technical Specification:

| | |
|--------------|--|
| Head | 5 pce Sliding Binocular Head, 30 ° Inclined, Interpupillary Distance 55-75mm |
| Eyepiece | WF10X |
| Objective | Achromatic 4X |
| | Achromatic 10X |
| | Achromatic 40X(S) |
| | Achromatic 100X(S, Oil) |
| Nosepiece | Quadruple Revolving Nosepiece |
| Stage | Double Later Mechanical Stage 140mm x 140mm, Travel range 75x50mm |
| Condenser | Abbe N.A.=1.25 with Iris Diaphragm & Filter |
| Focusing | Coaxial Coarse & Fine Focusing Knobs |
| Light Source | Base with Built-in Halogen Lamp 6V/20W, with Continue intensit Control, |
| | 110V/220V |
| Indicator | Demo. Indicator with Illuminator of 6V/3W |

Regulatory compliances



Corporate Social Responsibility



Analytical
Foundation

Analytical Foundation is a nonprofit organization (NGO) founded for the purpose of:

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various fields of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org
2. Improving quality of life by offering YOGA Training courses, Workshops / Seminars etc.
3. ANALYTICAL FOUNDATION aims to DETOXYFIFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Comps etc.

Reach us @



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide