

# PA - 3002

## Polarimeter Saccharimeter



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

## ►► Product Introduction

Polarimeters are used for the determination of the Optical Rotation of all Optically Active substances and used in testing Pharmaceuticals, Chemicals, Essential Oils, Sugars, Foodstuffs, Brewery products etc..

The Polarimeters and Saccharimeters are available in two series. The PA Series is suitable for observation tubes upto 220mm while PB series can accommodate tubes upto 400 mm. They offer an accuracy of 0.01° Angular Scale and 0.5° ISS.

An outstanding feature of Polarimeters is the incorporation of precision glass scale and scale plate with minutely finished rotating mechanism and Tropic proof highly sensitive optical system. The test result can be read precisely and directly without the aid of a vernier and its associated errors. The latest addition is use of an Integrated Lighting System, in place of the bulky and fragile Sodium Lamp, which gives the advantage of immediate illumination while saving table space and electrical power.

The advance optical system of the Polarimeters and Saccharimeters have been designed in order to minimize the destructive influences of a hot and humid climate and abrupt temperature changes on the sensitive optics and are suitable for use in Tropical conditions.

Superior quality Optics & highly sensitive Triple Shadow offers the following advantages

- Permanent focus with difficult and dark liquids
- High accuracy and repeatability along with strain free viewing
- Triple Shadow System rotates with the circle thus affording immediate indication of the direction of the material being examined.

## ►► SAMPLE TUBES

A wide selection of sample tubes, like Side Filling, Center Filing and Jacketed tubes are available. Sample tubes are available from 10mm/0.1dm, 25mm/0.25dm (both in stainless steel), 50mm/0.5dm , 100mm/1dm, 200mm/2dm, 220mm/2.2dm, 400mm/4dm. Jacketed tubes for precision temperature control of samples are available in 100mm 200mm lengths .Centre filling tubes are also supplied on request.

## ►► Technical Specifications

Instrument Name	Polarimeter	Polarimeter	Polarimeter	Polarimeter/ Saccharimeter	Polarimeter/ Saccharimeter
Model No.	PA-3002	PA-3001R (Standard)	PA-3001R (Deluxe)	PB- 3001R(Standard)/ PB-3004 (Standard)	PB- 3001R(Deluxe)/ PB-3004 (Deluxe)
Tube Length (Maximum)	220 mm	220 mm	220 mm	400 mm	400 mm
Shadow System	Half/Triple Shadow	Improved Triple Shadow	Imported Triple Shadow & High Definition Optical System	Improved Triple Shadow	Imported Triple Shadow & High Definition Optical System
Angular Scale	0-360 <sup>0</sup>	0-360 <sup>0</sup>	0-360 <sup>0</sup>	0-360 <sup>0</sup>	0-360 <sup>0</sup>
Angular Scale Accuracy	0.05 <sup>0</sup>	0.01 <sup>0</sup>	0.01 <sup>0</sup>	0.01 <sup>0</sup>	0.01 <sup>0</sup>
International Sugar Scale	±130 <sup>0</sup> Z	±130 <sup>0</sup> Z	±130 <sup>0</sup> Z	±130 <sup>0</sup> Z	±130 <sup>0</sup> Z
ISS Accuracy	0.1 <sup>0</sup> Z	0.05 <sup>0</sup> Z	0.05 <sup>0</sup> Z	0.05 <sup>0</sup> Z	0.05 <sup>0</sup> Z
Light Source	Sodium Lamp/ Integrated Light Source	Sodium Lamp/ Integrated Light Source	Sodium Lamp/ Integrated Light Source	Sodium Lamp/ Integrated Light Source	Sodium Lamp/ Integrated Light Source
<b>Standard Sample Tube</b> (Side Filling With Spare Window Glass and Washer)	200 mm: 01 No.	100mm : 01 No. 200mm : 01 No	100mm : 01 No. 200mm : 02 No	200mm : 01 No. 400mm : 01 No	100mm : 01 No. 200mm : 01 No 400mm : 01 No.
Size	27"X13 <sup>1/2</sup> "14 <sup>1/2</sup>	28"X13 <sup>1/2</sup> "14 <sup>1/2</sup>	28"X13 <sup>1/2</sup> "14 <sup>1/2</sup>	28"X13 <sup>1/2</sup> "14 <sup>1/2</sup>	28"X13 <sup>1/2</sup> "14 <sup>1/2</sup>
Weight	11 kg	12 kg	12 kg	15 kg	15 kg
Cold/Hot water circulating baths for sample tube temperature control available as additional accessory					

## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

**HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

**Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

**AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

**Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

**Instruments** :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## Our Products & Technologies



UV/VIS  
Spectro 2080+  
Double Beam



Infra FTIR



Optima Gas  
Chromatograph  
3007



Optima Gas  
Chromatograph  
2979 Plus



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Particle  
Counter



Optical Emission  
Spectrophotometer



DSC/TGA



Semi Auto Bio  
Chemistry Analyzer



HEMA 2062  
Hematology  
Analyzer



Micro Plate  
Reader/Washer



URINOVA 2800  
Urine Analyzer



Total Organic  
Carbon 3800



Fully Automated  
CLIA



NOVA-2100  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



TOC  
Analyzer



Laser Particle  
Size Analyzer



Ion Chromatograph



Water purification  
system



## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



**HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)**

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W: [www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)  
[www.ais-india.com](http://www.ais-india.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide