

VET - 3026A

Portable Color Doppler Ultrasonic System




EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Ergonomic Design




►► Features

- 2D, CFM, M, PW, CW, CMM, DPDI, Biopsy Function
- THI, SRI, SCI, TSI, FCI, EFOV, HR flow, B-Steer
- Support local and global amplification
- Support PW auto trace and calculation
- Preposed USB ports, supporting multi-data transmission
- Dual Live and multi-synchronization technologies

►► Accessories






Brilliant choice for those who require outstanding image quality, excellent portability and intelligent workflow.

►► **uSeed Platform**

M5 facilitates accurate diagnosis based on uSeed metadata beamforming Imaging technology platform. This brand-new platform is developed by CPU + GPU heterogeneous computing technology. The wave in uSeed platform are saved as metadata for postprocessing.




►► **Applications**

M5 VET has excellent performance and different applications according to user's demands: Felidae, Canine, Ovine, Equine, Bovine, Swine, User-defined
Dedicated software package: Comprehensive Measurement Package, Vet-specific Presets, Vet-specific Comments, Vet-specific Body Marks.

►► Technologies


Pulse Inverse Harmonic Imaging PIHI

PHPI reduces the distortion generated by the fundamental wave, and greatly improve the signal-to-noise ratio.



Speckle Reduction Imaging

SRI technology significantly reduce the unwanted speckle, provides enhance tissue boundry and tissue echo for confident clinical diagnosis.



SRI OFF

SRI ON

Spatial Compound Imaging

SCI provides better contrast resolution, reduces the speckle noise and smoothes the imaging of homogeneous tissue.



Traditional



Spatial Compound

Gallery



Canine-colon



Canine-color



Canine-Fetal heart



Canine-Fetal



Canine-Fetal-5W-Cord



Canine-gallbladder



Canine-uterus body



Felidae-B+C



Felidae-Bladder



Felidae-Descending colon



Canine-Left kidney



Canine-Liver& Gallbladder

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide