

VR - 18L/28L

Portable/Vehicle Refrigerator



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► **Portable lightweight design**

1. ABS Engineering Plastic Shell, 304 stainless steel lining, optional built-in wire basket, easy classification storage.
2. Standard Aluminum Alloy handle, strap fastening, shoulder back, hand can be easy to be fixed in the car trunk. (VR-18L)

►► **German technology, Quality of craftsmanship**

1. AC POWER SUPPLY VOLTAGE AC 220V ~ , Vehicle power supply voltage DC12/24V, perfect battery protection.
2. Germany SECOP DC compressor, Germany EBM Fan, micro channel condenser, high heat transfer efficiency, stable operation, especially suitable for vehicle transport environment.
3. The air-cooled circulating fan is arranged in the box, the temperature is stable and uniform, and the fluctuation is small, Built-in LED lighting, easy to operate.
4. High-pressure foaming insulation layer, Non-fluorine refrigerant.
5. Temperature setting range:- 20°C - 10°C, accuracy 0.1°C (ambient temperature 32°C).

►► High precision control system, more security guarantee

1. Embedded intelligent control system, 2.8-inch TFT touch screen, easy to operate, dynamic display running status.
2. Optional 4G wireless networking function, Docking cold-chain monitoring cloud platform, storage temperature monitoring anytime, anywhere, one more security protection.
3. Optional Bluetooth Printer, easy to print temperature history at any time.
4. Support USB Data Export Function, can export PDF, excel tables, support the generation of statistical reports or charts.

Model	External size (mm)	Capacity (L)	Temperature (°C)	Power (W)	(kg) Net Weight
VR-18L	500*350*350	18	10~ -20	42	10.5
VR-28L	500*350*440	28	10~ -20	50	12

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.


Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph




Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide