



Refrigerator (-20C)

▶▶▶ Technical Specification:

Model : 3450
Effective Capacity : 450 Liters
Temperature Range : -10°C ~ -25°C

▶▶▶ Specs:

- ▶ Vertical cabinet (upright), with two chamber rooms.
- ▶ Temperature range: from - 10°C to -25°C, can be set freely.
- ▶ Microprocessor temperature controller LCD display, to make clear the temperature information.
- ▶ ABS material interior, Outer body is painted steel board
- ▶ The safety lock on the door
- ▶ Retatable sure-grip door handle
- ▶ Perfect audible/visual alarm: High or low temperature alarm, system failure alarm
- ▶ 60 mm extra thick insulation
- ▶ The famous compressor (2unit), can control up and down temperature separately.
- ▶ Password protected configuration page
- ▶ 4 Casters provided for easy handling
- ▶ 12units drawers made of ABS material inside.
- ▶ Refrigerant as R134A, CFC free
- ▶ Power supply: 220V/50Hz 1 phase, can change as 220V 60Hz or 110V 50Hz or 110 V 60Hz
- ▶ Certificate: ISO9001, ISO14001, ISO13485

▶▶▶ Specifications:

Model	capacity	External size (W*D*H)mm	Inside size (like (W*D*H)mm	Input power	Weight (Nt/Gt)
3450	450liters	655*780*1845	530*626*1356	450W	120/140

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶▶ Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.