

RE-3002

Rotary Evaporator



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

RE- 3002

The rotary evaporators forms the ideal basis for efficient evaporation in your laboratory. Aspects such as safety, convenience and cost savings have been considerably taken into account during the development of the evaporators. This is reflected in various unique product advantages. Among others all evaporators are equipped with a separate ON/OFF switch for the heating function. It is also provided with easy lift, a simple touch operation to lift flask. With the rotary evaporators you can as well compose your individual modular concept, consisting of rotary evaporator, vacuum controller, vacuum pump and chiller. Distill your solvents reliably. It offers an intelligent solution to process your solvent evaporation easy and safely. Our vacuum distillation system makes evaporating easy for you and allows you to concentrate on your research.

Features:

- Comes with an integrated Combo-Clip
- Equipped with an LED rotation speed setting
- Comes with an electronic speed control rotation drive
- Has a flexible bath position
- Equipped with a vapor tube that can be easily removed for cleaning
- Comes with a safety protection switch to avoid overheating
- Manual lift type
- Sealing Ring (Vacuum Seal): Graphite filled PTFE
- Ceramic coated heating bath
- Equipped with a secondary protection cut off switch

RE- 3002A

This Rotary evaporators forms the ideal basis for efficient evaporation in your laboratory. They are suitable for standard distillation, for complex distilling processes as well as for automatic distillation. Aspects such as safety, convenience and cost savings have been considerably taken into account during the development of the evaporators. This is reflected in various unique product advantages. Among others all evaporators are equipped with a separate ON/OFF switch for the heating function. It is also provided with easy lift, a simple touch operation to lift flask. With the rotary evaporators you can as well compose your individual modular concept, consisting of rotary evaporator, vacuum controller, vacuum pump and chiller. Distill your solvents reliably. It offers an intelligent solution to process your solvent evaporation easy and safely. Our vacuum distillation system makes evaporating easy for you and allows you to concentrate on your research.

Features:

- Comes with an integrated Combo-Clip
- Equipped with an LED rotation speed setting
- Comes with an electronic speed control rotation drive
- Has a flexible bath position
- Equipped with a vapor tube that can be easily removed for cleaning
- Comes with a safety protection switch to avoid overheating
- Automatic lift type
- Sealing Ring (Vacuum Seal): Graphite filled PTFE
- Ceramic coated heating bath
- Equipped with a secondary protection cut off switch

RE-3005A

This rotary evaporators forms the ideal basis for efficient evaporation in your laboratory. They are suitable for standard distillation, for complex distilling processes as well as for automatic distillation. Aspects such as safety, convenience and cost savings have been considerably taken into account during the development of the evaporators. This is reflected in various unique product advantages. Among others all evaporators are equipped with a separate ON/OFF switch for the heating function. It is also provided with easy lift, a simple touch operation to lift flask. With the rotary evaporators you can as well compose your individual modular concept, consisting of rotary evaporator, vacuum controller, vacuum pump and chiller. Distill your solvents reliably. It offers an intelligent solution to process your solvent evaporation easy and safely. Our vacuum distillation system makes evaporating easy for you and allows you to concentrate on your research.

Features:

- Comes with an integrated Combo-Clip
- Equipped with an LED rotation speed setting
- Comes with an electronic speed control rotation drive
- Has a flexible bath position
- Equipped with a vapor tube that can be easily removed for cleaning
- Comes with a safety protection switch to avoid overheating
- Automatic lift type
- Sealing Ring (Vacuum Seal): Graphite filled PTFE
- SS316 heating bath
- Equipped with a secondary protection cut off switch

Technical Specifications:

Description	RE- 3002	RE- 3002	RE- 3002
Lift Type	Manual	Manual	Manual
Condenser Type	Vertical/Horizontal	Vertical/Horizontal	Vertical
Cooling Surface Area	1500 cm ²	1500 cm ²	2000 cm ²
Rotation Speed (RPM)	20-320	20-320	20-320
Vapour Tube	Single piece with Sleeve for Easy Removal	Single piece with Sleeve for Easy Removal	Single piece with Sleeve for Easy Removal
Integrated Combi-Clip	Yes	Yes	Yes
Suitable flask size (ml)	50-5000 ml	50-5000 ml	50-5000 ml
Rotation speed setting	LED Display	LED Display	LED Display
Motor Power (W)	60 (Brushless DC)	60 (Brushless DC)	60 (Brushless DC)
Heating Capacity (W)	1300	1300	1300
Total Power (W)	1360	1360	1360
Temperature Range (°C)	20 – 210	20 – 210	20 – 210
Temperature Accuracy(°C)	1	1	1
Overheat Protection	Safety cutoff over +5 (°C)	Safety cutoff over +5 (°C)	Safety cutoff over +5 (°C)
Temperature Setting	Digital LED	Digital LED	Digital LED
Temperature Control	Electronic PID Control	Electronic PID Control	Electronic PID Control
Secondary Protection Cut off	250°C	250°C	250°C
Material Heating Bath	Ceramic Coated	Ceramic Coated	SS316L
Volume Heating bath (L)	4.5L-Dia 255mm Suitable 50- 5000ml	4.5L-Dia 255mm Suitable 50- 5000ml	5L Rectangular
Bath Position	Flexible to extend, to suit flask size	Flexible to extend, to suit flask size	Flexible to extend, to suit flask size
Sealing ring	Graphite filled PTFE	Graphite filled PTFE	Graphite filled PTFE
Power Supply	230 VAC, 50Hz	230 VAC, 50Hz	230 VAC, 50Hz

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



LCMS



Automated Prep-Flash Chromatography system



Maldi TOF



Optima Gas Chromatograph



Flash Chromatograph



DAC Column



GCMS 3068




UHPLC



HPTLC



Ion Chromatograph



Production HPLC



Helium Mass Spectrometer Leak Detector



Column



DLS



Water purification system

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

ANALYTICAL FOUNDATION is a Nonprofit Organization (NGO) for the purpose of:



Analytical Foundation


1. Research & Innovation **Scientist's awards & International Conference visit sponsorship / QC Professional Awards:** Quality life is possible by innovation only and the innovation is possible by research only, hence **ANALYTICAL FOUNDATION** is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research /testing/ publication at info@analyticalfoundation.org

2. Improving quality of life by offering **YOGA Training courses, Work shops / Seminars** etc.

3. **ANALYTICAL FOUNDATION** aims to DETOXIFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Naturecure, Health Care, Ara, Events, Camps etc.

▶▶ A portion of profit of **ANALYTICAL GROUP** Company's goes into **ANALYTICAL FOUNDATION** for the aboce noble cause.

▶▶▶ Reach us @



Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice

HPLC Solutions	MultipleLabs	Infinite Multiplelabs LLP	Analytical Bio-Med	Analytical Foundation (Trust)
----------------	--------------	---------------------------	--------------------	-------------------------------

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E :- info@analyticalgroup.net
info@hplctechnologies.com
info@multiplelabs.com

W:- www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide
USA| UK| Australia| Dubai| Malaysia