

SDP-3002

Statistical Data Printer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► ATL STATISTICAL PRINTER is ideal for routine laboratory application with following features.

ISO/GLP compliant recording capability for printouts.

GLP : Good Laboratory Practice Password-protected menu lock.

We have the parameter pertaining to the ambient weighing conditions printed before (GLP header) and after (GLP footer) the values of a weighing series. These parameters include:

GLP header :

- Date
- Time at beginning measurement
- Balance manufacturer
- Balance model
- Balance serial number
- Software version number
- Identification number of the current sampling operation

GLP Footer:

- Date
- Time at end of measurement
- Field for operator signature
- ID designations are configured as follows setup: Printout : Identification codes.
- ISO/GLP compliant recording capability for printouts.

01:	50.000 g
02:	100.000 g
03:	150.000 g
04:	200.001 g
<hr/>	
n	4
\bar{x}	125.0003 g
s	64.5501 g
srel	51.63 %
Σx	500.001 g
min	50.000 g
max	200.001 g
Diff	150.001 g
<hr/>	
25-Mar-2010	3:33:46

22-Mar-17	11:13AM
<hr/>	
Calibration:	Internal
W-ID	*****
Temperature	29.378°C
Set	+ 200.0000g
Diff.	-0.0031g
<hr/>	
Internal Cal Done	
Diff.	0.0000g
22-Mar-17	11:13AM
Name :	*****



Model	SDP-3002
Features	Date & Time (Password Protection), Statistic & totalization function
Printer	Dot Matrix Printer
No.of Characters per line	24 Character
Printing Speed	Faster than 1 line per second
Paper Type	Normal Paper Roll, Wood Free Paper
Paper Size	Width:57.5 ± 0.5mm.
Power Supply	AC Adapter, 100-270V AC, 50Hz-60Hz
Interface	RS 232C
Bound Rate(Selectable)	1200,2400,4800,9600,19200
Parity(Selectable)	None, Space, Odd, Even
Data Bit	8 Bit
Stop Bit	1
Working Temperature	5°C to 50°C, relative humidity : 10%-80% (no condensation)
Dimensions(mm/inch)(L x W x H)	210 x 135 x 88.5/ 8.30" x 5.30" x 3.50
Net Weight/Gross weight	2.00 kgs / 2.5 kgs
Packing Size(mm/inch)	270 x 200 x 175/10.62" x 7.8" x 6.90"

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.