

Suction Machine

SP-3001



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Specifications

- High vacuum suction unit, run on electricity with two suction jars of 3 liters capacity each (Suction capacity 50-60L/min).
- Maximum vacuum 650mm Hg
- Auto cut off device for preventing entry of fluid in pump.
- 4 castors
- Fast and efficient jar change facility.
- Easy access and controls.
- It be heavy duty and noiseless.
- Be able to create desired maximum vacuum in least possible time.
- Light and manoeuvrable.
- One plastic suction jar cover, steam sterilizeable to be provided extra.
- Two extra suction jar (plastic) of capacity 3 liters be provided.
- Be quoted along with accessories like lid, tubing etc with the equipment to make the unit functional.
- Mechanical overflow protection system
- Max noise level > 50 db
- Power input 220-250VAC, 50Hz

►► Regulatory compliances



►► Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

►► Reach us @



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents /technical specifications of the catalogue without prior notice.