

TOC detection of purified pharmaceutical water and Water For Injection(WFI)

TOC 3600i



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES


Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net


► Test Methods and Samples Management

Meet the GMP and Computerized System Validation requirements




Electric Record and Signature

- 3-levels password authorization
- DI1500 will establish complete "Users List"
- Every user has the unique and uncreateable user name
- Help the user prepare the SOP



Data Audit Trail

- All operations are kept in data log
- Data log can be looked up
- Data log cannot be modified and deleted



Raw Data Tracibility

- All data is stored in real time
- Follow the data tracibility regulations
- Raw data backup and recovery function
- All data can be read after software updated

Online Modular IP65 Protection



►► Features

- Complying with the requirements of 21CFR Part 11 Powerful software to improve efficiency
- Fully automatic operation
- UV lamp oxidation, no need catalyst and carrier gas
- Autosampler for option
- Online modular for option
- IQ,OQ validation documents are available

►► Parameters

Power	100-240V 50/60Hz	Power Consumption	100W
Readability	0.001mg/L	Detection Accuracy	±3%
Detection Range	0.001mg/L-1.500mg/L	Water Conductivity Required	0.55-2.00 us/cm
Response Time	less than 15 min	Analysis Time	4 min.
Ambient Temperature	10-60°C	Samples Temperature	1-95°C



Autosampler AS-3670



Autosampler AS-3618



►►► Regulatory compliances



►►► Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

4. buy our products and be a part of social cost of analytical foundation.

►►► Reach us @



Analytical[®] Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.