

Temperature water Bath with digital control



Constant temperature water bath is equipped with a microprocessor intelligent control system. With circulating pump, it can let the heated liquid flow off from the tank and thus establish a second constant-temperature field. Available is a large-screen LCD display. It is quick to set the temperature by touching the soft keys and easy to operate. The Micro-Computer can amend deviation of temperature measurements.

Features:

➤ Temperature Range °C	Ambient+ 5-100°C
➤ Temperature Stability °C	±0.1°C
➤ Temperature Accuracy °C	±0.1°C
➤ Display Resolution (mm3)	0.1/0.01
➤ Inner Tank size (mm3)	280× 190 ×200
➤ Depth (mm)	200
➤ Top Opening Dimension (mm)	180×160
➤ NW (kg)	14

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide