

# WATER DISTILLATION SYSTEM

(Fully Automatic Cabinet Model)

## SPECIAL FEATURES

Its a fully automatic water distillation system with advanced controlled technology and can be used in R&D laboratories,hospitals etc.

It is manufactured with optional quartz glass or borosilicate glass.

Latest unique Microprocessor technology with **Touch Screen TFT** or **Graphical LCD Display** makes it elegant in look and user friendly in operation.

Borosilicate condenser and boiler is fit with quartz heater which is removable and screw threaded connector to clean the system when scaling occurs.

Cooling water recirculation for energy saving

Automatic Detection of Water Flow

Level sensor to prevent Heater from Damaging

Vessel Full Level Sensor

Timer Based Opeartion (Double Stage Model)

Audio / Visual Alarm if no water supply

## DOUBLE STAGE

## SINGLE STAGE

### Technical Specification

Model	WDS-30-S2	WDS-30 -D2	WDS-30-S4	WDS-30 -D4
<b>Output Capacity</b>	2	2	4	4
<b>STAGE LEVEL</b>	SINGLE	DOUBLE	SINGLE	DOUBLE
<b>Cooling water</b>	2 Ltr/Min	1.5 Ltr/Min	2-3 Ltr/Min	2-4 Ltr/Min
<b>Power Consumption</b>	1.5KW	1.5 *2 KW	3.0 KW	3.0*2KW



## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



**Analytical**  
**Foundation**

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



 **Analytical**<sup>®</sup>  
Technologies Limited

HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.ais-india.com](http://www.ais-india.com)  
[www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide