

Analytical Technologies

XRF 3005



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Details

Based on the mature technology of X-ray fluorescence analyzer and the latest technological development, and on the performance of the existing international commercial handheld X-ray fluorescence analyzer, adding new elements into it on the basis of the first, second, and third generation, Analytical Technologies has developed the intellectual portable fourth generational handheld analyzer XRF series, which can satisfy the market requirements and possess independent intellectual property rights.

The new generational XRF series are designed according to the applications of on-site X-ray analysis in the wild, featuring small and light body, which can be held in hands by common people when measuring. It is super small, super light, super beautiful, super safe, super convenient, super long standby time, super waterproof, super precise and super fast. Digital multi-channel technology make the instrument have better detection limit, better stability, and wider application fields.

Products include RoHS analyzer 3003XRF, alloy analyzer 3005XRF, mineral analyzer 3007XRF, and heavy metals in soil analyzer 3009XRF, which can be applied for in-situ material testing, for on-site on-line analysis of various elements during the process from manufacturing control to finished product testing, for quick identification of alloy number and reliable authentication, for on-site heavy metals testing in soil, and for various kinds of geological and ore samples on-site analysis of multi-elements.

This series of products has passed strict testing and inspection, and all indexes are accorded with the related technological requirements, reaching the international advanced level. Products have been sold both home and abroad in batch, getting favourable comments from customers, breaking the monopoly of foreign manufacturers. Performance Advantage Perfect performance as Introducing three core technology of small power window integrated miniature X-ray tube, large dimension beryllium window electric-cooling SDD detector (the best detector in the world), and miniature digital signal multi-channel processor, greatly reduce the testing time and test deviation, and improve the detection precision.

Small body, easy to carry and operate. Small body. Easy to carry. Convenient for wild work. Can operate on-site and in-situ analysis at anytime or anywhere.

Nondestructive detection. No destruction to samples.

Fast detection. Instrument not only can operate fast detection with handheld type within 1-2 seconds, but also can operate long time precise detection with type, having the similar precision as the laboratory detection after 10 seconds.



Newly added light elements detection function Compared to handheld X-ray fluorescence analyzer of the third generation, it is added the gas charging system, which can charge helium at ordinary pressure. Thus it can detect the elements from Mg, greatly expanded the detecting range of elements, satisfying the requirements of customers according to light elements detection.

Direct testing It can detect the analytes directly, no need of samples. It can analyse various kinds of samples, including electronic products, alloy samples, geological and mining, soil, rock, residues, small solid particles, liquid sediments etc.

HD camera, more accurate testing Built-in HD camera, can observe the testing position at any time, which is very important of mineral sample testing.

Easy for calibration of deviation Multi-testing modes installation and free adding of infinite modes, coordinating with automatic testing mode matching function, it achieves the easy test for one key. Built-in intensify calibration methods, it can calibrate the deviation caused by different geometry shapes and inhomogeneous structure density.

Professional software, easy to operate RoHS elements analyzer, alloy analyzer, mineral analyzer, and heavy metals in soil analyzer all are installed with professional software. Brand new software interface and core, combined FP with EC, it has wider application fields.

The collimator filter, improving 3 times of efficiency On the basis of the combination of collimator and filter, increase them from 4 groups to 12 groups, increasing 3 times.

Faster data transmission Embedded Window CE system, HD touch screen (resolution 640*480), digital multi-channel technology, and SPI data transmission, effectively improve the data transmission ability and counting ability, mastering testing data in every environment.

Radiation protection, care for health Triple safety protection function, automatic sensation, automatically shut-down of X light pipe within 2 seconds with no sample in testing; The radiation level is far lower than the international safety standard when working, with no radiation leakage; accompanied with test safety cover.

Strong power, convenient to charge Lithium battery, with maximum capacity of 7800 mAH, can continuously work for 8 hours, whose endurance are 2 times longer compared to last three generations. And equip with wide voltage ac charger or vehicular charger, which ensure the test at any place or any time.

Warning instructions prevent for mistaken operation Equipped with warning instruction system, green light flashes when the power is on and the yellow light flashed when testing, which prevent for mistaken operation.



Multiple protection and durableThe instrument has waterproof and dustproof function and can continuously work under high temperature and humidity. Its protection boxes are manufactured with high strength military material ,which enjoy with well moistureproof, shockproof, and pressure resistant functions

Performance index

1. Model:The 4th generation EDX analyzer XRF series.
2. Analysis method:Energy dispersive X-ray fluorescence analysis method.
3. Measuring range of elements: Mg to U.
4. Simultaneously detect elements:Simultaneously detect tens of elements.
5. Processor and RAM:CUP: 667MHz, RAM:256M , Maximum expanded storage: 32G, Standard configuration: 2G , can store large amounts of data.
6. Content range: ppm~99.99%.
7. Testing time :3-30 seconds.
8. GPS 、WIFI:Built-in system.
9. Battery time: Lithium battery,which can be charged, with maximum capacity of 7800mAH, can continuously work for 8 hours; Equiped with wide voltage(110V-220V) general adapter, can work under alternating current.
10. Testing object:Solid, liquid , powder.
11. Detector:25mm² 0.3mil, SDD detector.
12. Detector resolution: Lowest resolution can be 139eV .
13. Excitation source:40KV/100uA-Ag anode window miniature X light tube and high voltage source.
14. Collimator and filter:Collimator of 4.0 or 2.0 diameter, automatic switch of 6 types filter groups. 12 kinds of groups, world's most compound mode, can satisfy various kinds of samples testing.
15. Video system:HD camera.
16. Screen :TFT-LCD touch screen, resolution 640*480.
17. Detection limit:Lowest detection limit accounts to ppm level.
18. Testing window:12mm.20.Safety :Self-contained password manager mode.
19. Gas charging system: Helium charging at ordinary pressure system.
20. Data transmission: Digital multi-channel technology, SPI data transmission, quick analysis, high counting rate, waterproof miniature USB,which can be connected to computer.
21. Humidity :≤90%.
22. Temperature: -20°C ~ +50°C.
23. Size :234×306×82mm (L×H×W) .
24. Weight:1.9Kg (with battery) , 1.6Kg (without battery).

3003 XRF

The instrument is widely used in toy safety testing, electrical and electronic and RoHS testing etc. .



HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing : HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings : We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC : AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations : Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments : Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies

						
UV/VIS Spectro 2060+ Double Beam	FTIR	Gas Chromatograph 3000	Gas Chromatograph 2979 Plus	Flash Chromatograph	Atomic Absorption Spectrophotometer	Liquid Partical Counter
						
Optical Emission Spectrophotometer	DSC/TGA	NOVA 2020 plus Automated Bio Chemistry Analyzer	HEMA 2020 Hematology Analyzer	Micro Plate Reader/Washer	URINOVA 2800 Urine Analyzer	Total Organic Carbon
						
Fully Automated CLIA	NOVA Basic Semi -Auto Chemistry Analyzer	PCR/Gradient PCR/ RTPCR	Blood Gas Analyzer	Random access Analyzer for immunoassay Proteins & clinical chemistry	Semen Analyzer	Water purification system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.