

DISINTEGRATION TESTER 6 3000 Series



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Features :

- Complies with USP, IP and Ph. Eur. specifications
- Supports 6 USP baskets
- Each pair of baskets controlled by independent timer and motor
- Loading and a swivel free basket movement
- Green illumination for better visibility of disintegration process
- Programmable temperature controller
- Magnetically coupled pump for water circulation in the bath
- Digital Indicator for Temperature with Optically isolated Electronic heater driver with water bath
- Smooth Maintenance free drive
- Programmable Timer in mins/secs
- Sleek & Rugged..
- Soft membrane Key pad & Splash proof Panel
- Ergonomically designed Control Panel
- Gearless Drive
- Swivel free Basket Movement
- Green light illumination increases visibility

►► Specifications :

No. of stroke	30 +0.5 stroke/minute
Stroke height	55 mm ±2 mm
Motor	3 (DC stepper motor of 6V DC and 3.5 kg cm Torque)
No. of baskets	6 nos. (Each having 6 test tubes)
Water bath	6 mm thick transparent acrylic
Heater	400 watts, S.S. 304
Temperature range	20°C to 39.9°C
Temperature resolution	0.1°C
Temperature control accuracy	±0.5 °C
Timer	Range 1- 1 sec. to 99 min. & 59 sec.
	Range 2- 1 min. to 9 hrs & 59 min.
Timer resolution	Range 1- 1 sec.
	Range 2- 1 min
Power supply	220/230 VAC, 50/60 Hz, 500 VA
Display	4 Digit LED each
Dimensions (L x W x H)	350 x 390 x 510 mm (approx.)
Weight	Net 23 kgs and Gross 28 kgs

►► HPLC Servicing, Validation, Trainings and Preventive Maintenance :

- HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
- Trainings :We also take up preventive maintenance to reduce downtime of HPLC's Trainings.
- AMC's/CMC : We offer user training both in-House and at customer sites on PLC principles, operations, troubleshooting.
- Validations : We have protocols for carrying out periodic Validations as per GLP/ GMP/U SFDA norms.
- Instruments : We offer instruments/Renting Services Modules like pumps, detector etc. on Rent.



►► About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufaturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

►► Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

►► Regulatory compliances



►► Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

►► Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents /technical specifications of the catalogue without prior notice.