

Multi - Parameter Monitor ACMS 8600



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Description

ECG

Lead type	: 5-lead
Input	: RA;La;RI;LL;V
Sweep Speed	: 12.5mm/s, 25mm/s, 50mm/s
Accuracy	: ± 1 bpm or $\pm 1\%$ which is grater
Protection	: with Stand 4000VAC/50Hz Voltage in Isolation against electrosurgical and defibrillation
S-T Detection	: Yes, Arrhythmia analysis: Yes
Alarm	: Yea, audible and visual alarm, alarm events recallble

Standard Configuration

- ECG, NIBP, SPO2
- 8.4" color TFT Display
- Suitable for adult, Pediatric and neonatal patient
- Basic Parameters (ECG, SPO2, NIBP) in a durable case for bedside monitoring and transport monitoring
- Audible and visual alarms with adjustable alarm ranges
- Networkable CONTECT central monitoring system
- Powerful data management and storage capacity

NIBP

Method	: Oscillometry
Operation Model	: Manual/ Automatic/ SATA
Measurement Unit	: mmHg/KPa selectable
Measurement type	: systolic, Diastolic, And Mean
Over- Pressure Protection	: Yes

SPO2

Measurement Range	: 0-100%
Accuracy	: \pm digit (70-100%); 0-69% Unspecified.


Optional

RESP, 2-TEMP, 2-IBP, ETCO2, Network, Thermal Printer, Pace Maker

ETCO2 MODULE (Side stream)



IBP Module



ETCO2 Sampling Cube



IBP Transducer



Printer



IBP Cable



Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.