

Pipette Equipment PE - 3043



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Specifications

Used for determining the particle size distribution of very fine soils.

Basic items to perform the test are the pipette which is held by the Pipette stand for raising or lowering precisely the pipette to the required level without disturbing the suspension under test. The sedimentation cylinder is submerged in the Constant temperature water bath.

- Dimensions: 600x300x380 mm
- Weight approx.: 10.5 kg

▶▶ **HPLC Servicing, Validation, Trainings and Preventive Maintenance :**

- HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
- Trainings :We also take up preventive maintenance to reduce downtime of HPLC's Trainings.
- AMC's/CMC : We offer user training both in-House and at customer sites on PLC principles, operations, troubleshooting.
- Validations : We have protocols for carrying out periodic Validations as per GLP/ GMP/U SFDA norms.
- Instruments : We offer instruments/Renting Services Modules like pumps, detector etc. on Rent.



▶▶ **About Analytical Technologies**

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

▶▶ **Our Products & Technologies**

						
UV/VIS Spectro 2060+ Double Beam	FTIR	Gas Chromatograph 3000	Gas Chromatograph 2979 Plus	Flash Chromatograph	Atomic Absorption Spectrophotometer	Liquid Partical Counter
						
Optical Emission Spectrophotometer	DSC/TGA	NOVA 2020 plus Automated Bio Chemistry Analyzer	HEMA 2020 Hematology Analyzer	Micro Plate Reader/Washer	URINOVA 2800 Urine Analyzer	Total Organic Carbon
						
Fully Automated CLIA	NOVA Basic Semi-Auto Chemistry Analyzer	PCR/Gradient PCR/ RTPCR	Blood Gas Analyzer	Random access Analyzer for immunoassay Proteins & clinical chemistry	Semen Analyzer	Water purification system

▶▶ **Regulatory compliances**



▶▶ **Corporate Social Responsibility**

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶ **Reach us @**



HPLC Solutions	MultipleLabs	Analytical Bio-Med	Analytical Distributors	Analytical Foundation (Trust)
Corporate & Regd. Office: Analytical House, # E67 & E68, Ravi Park, Vasna Road, Baroda, Gujarat 390 015. INDIA	T: +91 265 2253620 +91 265 2252839 +91 265 2252370 F: +91 265 2254395	E: info@hplctechnologies.com info@multiplelabs.com info@analyticalgroup.net info@analyticalbiomed.com	W. www.ais-india.com www.analycalgroup.net www.hplctechnologies.com www.multiplelabs.com	Sales & Support Offices: across the country : Distributors & Channel partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.