Multiplelabs™
Turnkey Laboratories Solutions

Analytical®
One Brand, Many Solutions

ICF-3307

Ice Flaking Machine



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

ICF-3307

Features

- Hygienic flake ice production
- Compact space-saving design
- Automatic operation
- Stainless-steel construction
- Low power consumption
- Built-in safety protection
- Continuous ice production
- Easy maintenance
- Ideal for laboratory and industrial use

Model	ICF-3307
Capacity of Ice Flaker	300kg/24h
Storage Capacity	120kg
Storage Bin	ABS Plastic
Cooling Method	Air (inside R290A refrigerant)
Ambient Temperature	Below 25°C for optimum working
Construction	Stainless Steel
Noise Level	≤52 dB
Water Supply	3/4" BSP
Refrigerant	R290A
Power Supply	AC 220V/50 Hz
Input Power	1360 W
Dimensions (mm)	560x870x1580mm
Gross Weight	140 kg

Analytical®
Technologies Limited

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly the most affordable cost.

Trainings : We also take up preventive Maintenace to reduce downtime of HPLC's Trainings.

AMC's/CMC : We offer user training both in-House and at customer sites on HPLC principles, operations, troubleshooting.

Validations : We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments : We offer instruments/Renting Services Modules like pumps, detector, etc. on Rent.





About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufaturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



LCMS



Automated Prep-Flash Chromatography system



Maldi TOF



Optima Gas Chromatograph



Flash Chromatograph



DAC Column



GCMS 3068



UHPLC



HPTLC



Ion Chromatograph



Production HPLC



HPLC



Column



DLS



Water purification system

Regulatory compliances

Corporate Social Responsibility

Analytical Foundation

ANALYTICAL FOUNDATION is a Nonprofit Organization (NGO) for the purpose of:

1. Research & Innovation **Scientist's awards & International Conferance visit sponsership / QC Professional Awards**: Quality life is possible by innovation only and the innovation is possible by research only, hence **ANALYTICAL FOUNDATION** is committed to identify such personallities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research /testing/ publication at info@analyticalfoundation.org

2. Improving quality of life by offering **YOGA Training courses, Work shops / Seminars** etc.

3.**ANALYTICAL FOUNDATION** aims to DETOXIFY human minds, souls and body by means of Yoga, Meditation, Ayurveda,Naturecure, Health Care, Ara, Events, Camps etc.

A portion of profit of **ANALYTICAL GROUP** Company's goes into **ANALYTICAL FOUNDATION** for the aboce noble cause.

Reach us @



Analytical® Technologies Limited

Subsidiaries & Partners

USA Africa Russia Australia UK Malaysia China Dubai

SCAN ME

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370

E :- info@analyticalgroup.net
info@hplctechnologies.com
info@multiplelabs.com
analyticalfoundation@gmail.com
analyticalindia@gmail.com

W:- www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.analyticalfoundation.org

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide
USA| UK| Australia| Dubai| Malaysia

Company has copyrights to change the infromation of Brochure/Catalouge.