Multiplelabs™
Turnkey Laboratories Solutions

Analytical®
One Brand.... Many Solutions.

RC-3430

Refrigerated Micro Centrifuge



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

RC-3430

Features

- Bench top refrigerated micro centrifuge with max speed 20000 RPM.
- Heavy duty C.R.C. construction finished with powder coating
- Stainless Steel centrifugal bowl for extra protection and easy cleaning.
- Fitted with well balanced brushless motor with electronic drive.
- LCD Display for set time, run time, set and run rpm, rcf
- Digital Timer 1-99 min., lid lock and fast spin option and rotor Identification.
- Dynamic brake (automatic) and imbalance detector are provided.
- Selection of 3 acceleration + 3 deceleration profile.
- Last test parameter recall and displays automatically on LCD
- Microprocessor based. Programmable Speed Regulator from 2000 to 20000 in interval of 100 rpm
- P.I.D. Electronic digital temperature controller cum indicator is fitted to control the centrifuge temp. from room temperature upto -20°C.
- Supplied with complete , instruction manual, cord & plug
- Works on 230±10V/50Hz power supply.

Rotor Type	Rotor Capacity	Rotor Speed (RPM)
Angle rotor	12 x 1.5/2.2 ml	20000
Angle rotor	24 x 1.5/2.2 ml.	16000
Angle rotor	8 x 15 ml.	13000
Angle rotor	6 x 50 ml.	13000
Angle rotor	4 x 100 ml	13000
Angle rotor	4 PCR STRIPS OF 8 X 0.2ml	13000

Analytical®
Technologies Limited

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly the most affordable cost.

Trainings : We also take up preventive Maintenace to reduce downtime of HPLC's Trainings.

AMC's/CMC : We offer user training both in-House and at customer sites on HPLC principles, operations, troubleshooting.

Validations : We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments : We offer instruments/Renting Services Modules like pumps, detector, etc. on Rent.





About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufaturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies


LCMS


Automated Prep-Flash Chromatography system


Maldi TOF


Optima Gas Chromatograph


Flash Chromatograph



DAC Column


GCMS 3068


UHPLC


HPTLC


Ion Chromatograph


Production HPLC


HPLC


Column


DLS


Water purification system

Regulatory compliances

Corporate Social Responsibility

Analytical Foundation

ANALYTICAL FOUNDATION is a Nonprofit Organization (NGO) for the purpose of:

1. Research & Innovation **Scientist's awards & International Conferance visit sponsership / QC Professional Awards**: Quality life is possible by innovation only and the innovation is possible by research only, hence **ANALYTICAL FOUNDATION** is committed to identify such personallities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research /testing/ publication at info@analyticalfoundation.org

2. Improving quality of life by offering **YOGA Training courses, Work shops / Seminars** etc.

3.**ANALYTICAL FOUNDATION** aims to DETOXIFY human minds, souls and body by means of Yoga, Meditation, Ayurveda,Naturecure, Health Care, Ara, Events, Camps etc.

A portion of profit of **ANALYTICAL GROUP** Company's goes into **ANALYTICAL FOUNDATION** for the aboce noble cause.

Reach us @



Analytical® Technologies Limited

Subsidiaries & Partners

USA Africa Russia Australia UK Malaysia China Dubai

SCAN ME

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370

E :- info@analyticalgroup.net
info@hplctechnologies.com
info@multiplelabs.com
analyticalfoundation@gmail.com
analyticalindia@gmail.com

W:- www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.analyticalfoundation.org

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide
USA| UK| Australia| Dubai| Malaysia

Company has copyrights to change the infromation of Brochure/Catalouge.